

Product Datasheet

**GSK3 Alpha + Beta (phospho-Tyr279+Tyr216)
antibody [Out of stock] (orb526446)**

Description

Rabbit polyclonal antibody to

Tyr279+Tyr216) conjugated to PE

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

PE

Tested

FC, IF

Applications
Immunogen

KLH conjugated Synthesised phosphopeptide derived from human GSK-3 Beta around the phosphorylation site of Tyr216 VS(p-Y)IC

Target

GSK3B

Preservatives

0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.

Form/Appearance

Liquid

Concentration

- 100 µg (in 100 µl): 1 mg/ml

Storage

Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Uniprot ID

P49841

Dilution Range

Flow-Cyt: 1:50-200, IF: 1:50-200

Expiration Date

12 months from date of receipt.