



Product Datasheet

CD133 antibody (orb526485)


Description	CD133 antibody
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	KLH conjugated synthetic peptide derived from human CD133 (508-552/865 aa)
Target	CD133
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	95 kDa
Uniprot ID	O43490
Dilution Range	WB=1:500-2000, Flow-Cyt=1µg/Test
Expiration Date	12 months from date of receipt.



The blue histogram is unstained cells(He...



Flow cytometry of Mouse kidney.Primary A...



Western blot:SP2/0 Cell;(Mouse) Colon ca...