



Product Datasheet

CD14 antibody (orb526496)


Description	CD14 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	KLH conjugated synthetic peptide derived from mouse CD14 (241-300/366 aa)
Target	CD14
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	35/40 kDa
Uniprot ID	P08571
Dilution Range	WB=1:500-2000, Flow-Cyt=1ug/Test
Expiration Date	12 months from date of receipt.



Flow
cytometric
analysis of
Mouse
thymus...



Western blot
analysis of
Lung (Mouse)
Ly...



Flow
cytometric
analysis of
Raw264.7
usi...