

Product Datasheet

CTNNA1 (phospho-Ser641) antibody (orb526521)

Description	CTNNA1 (phospho-Ser641) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human CTNNA1 around the phosphorylation site of Ser641 (DD(p-S)DF)
Target	phospho-CTNNA1 (Ser641)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	100 kDa
Uniprot ID	P35221
Dilution Range	IHC-P=1:100-500
Expiration Date	12 months from date of receipt.