

Product Datasheet

MAP3K7IP1 (phospho-Ser423) antibody (orb526523)

Description	MAP3K7IP1 (phospho-Ser423) antibody
Species/Host	Rabbit
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human MAP3K7IP1 around the phosphorylation site of Ser423 (AH(p-S)AS)
Target	phospho-MAP3K7IP1 (Ser423)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	56 kDa
Uniprot ID	Q15750
Dilution Range	IHC-P=1:100-500, WB=1:500-2000, IHC-F=1:100-500, IF=1:100-500
Expiration Date	12 months from date of receipt.