

Product Datasheet

CD45 antibody (orb526577)

Description	CD45 antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	Recombinant human CD45
Target	human CD45
Preservatives	0.01M TBS(pH=7.4)
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	Mouse IgG2a, k
Clonality	Monoclonal
MW	143 kDa
Uniprot ID	P08575
Dilution Range	Flow-Cyt=1ul/Test
Expiration Date	12 months from date of receipt.