


Product Datasheet

NGAL antibody (orb526578)

Description	NGAL antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant human NGAL
Target	NGAL/Lipocalin 2
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	20 kDa
Uniprot ID	P80188
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.



Western
blot:A549(Human)
Cell;Primary:An...