


Product Datasheet

NMUR1 antibody (orb526582)

Description	NMUR1 antibody
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from mouse NMUR1/GPR66 (1-80/415 aa)
Target	NMUR1/GPR66
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	47 kDa
Uniprot ID	Q9HB89
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.



Western blot:Testis
(Mouse);Primary:Anti...