



## Product Datasheet

### HRH1 antibody (orb52659)


|                            |                                                               |
|----------------------------|---------------------------------------------------------------|
| <b>Description</b>         | HRH1 antibody                                                 |
| <b>Species/Host</b>        | Rabbit                                                        |
| <b>Reactivity</b>          | Human                                                         |
| <b>Conjugation</b>         | Unconjugated                                                  |
| <b>Tested Applications</b> | ELISA, IF, IHC, WB                                            |
| <b>Immunogen</b>           | Recombinant Human Histamine H1 receptor protein (211-416AA)   |
| <b>Target</b>              | HRH1                                                          |
| <b>Preservatives</b>       | Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol   |
| <b>Form/Appearance</b>     | Liquid                                                        |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| <b>Note</b>                | For research use only                                         |
| <b>Isotype</b>             | IgG                                                           |
| <b>Clonality</b>           | Polyclonal                                                    |
| <b>Uniprot ID</b>          | <b>P35367</b>                                                 |
| <b>Dilution Range</b>      | WB:1:200-5000, IHC-P:1:20-200, IF:1:50-200                    |
| <b>Expiration Date</b>     | 12 months from date of receipt.                               |



Immunohistochemistry of paraffin-embedde...



Immunofluorescent analysis of HepG2 cell...



Immunohistochemistry of paraffin-embedde...