

Product Datasheet

FBXO18 antibody (orb526609)

Description	FBXO18 antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	KLH conjugated synthetic peptide derived from human FBXO18 (1-100/1043 aa)
Target	FBH1
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/1ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	118 kDa
Uniprot ID	Q8NFZ0
Expiration Date	12 months from date of receipt.