

Product Datasheet

CD34 antibody (orb526618)

Description	CD34 antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human CD34
Target	human CD34
Preservatives	0.01M TBS(pH=7.4)
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
MW	39 kDa
Uniprot ID	P28906
Dilution Range	ICC=1:100
Expiration Date	12 months from date of receipt.