



Product Datasheet

STAT6 antibody (orb526632)

Description	STAT6 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	Recombinant human STAT6 protein, around 700-800aa
Target	STAT6
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Recombinant
MW	94 kDa
Uniprot ID	P42226
Dilution Range	WB=1:500-1000, Flow-Cyt=2ug/Test, IHC-P=1:100-500, IHC-F=1:20-200, ICC=1:20-200, IF=1:20-200
Expiration Date	12 months from date of receipt.



Flow cytometry of
MCF7.Primary
Antibody ...



Western
blot:NIH/3T3(Mouse)
Cell;HeLa(Hu...