


Product Datasheet

PINK1 antibody (orb526643)

Description	PINK1 antibody
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC-Fr, IHC-P, WB
Immunogen	Recombinant human PINK1 protein
Target	PINK1
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
MW	63 kDa
Uniprot ID	Q9BXM7
Expiration Date	12 months from date of receipt.



Western blot: All lanes
: Anti-Pink1(115-...