



Product Datasheet

Akt1 (phospho-Ser473) antibody (orb526646)


Description	Akt1 (phospho-Ser473) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	FC, ICC, IHC-P, WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human Akt1 around the phosphorylation site of Ser473 (QF(p-S)YS)
Target	Phospho-Akt1 (Ser473)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Recombinant
MW	56 kDa
Uniprot ID	P31749
Dilution Range	WB=1:500-2000, IHC-P=1:50-200, ICC=1:100, Flow-Cyt=2ug/Test, IHC-F=1:50-200, IF=1:50-200
Expiration Date	12 months from date of receipt.



Flow cytometry of U-937.Primary Antibody...



ICC analysis of NIH/3T3 cell using Phosp...



IHC-P analysis of mouse lung using Phosp...