

Product Datasheet

BCL2L1 antibody (orb526657)

Description	BCL2L1 antibody
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	KLH conjugated synthetic peptide derived from human BCL2L1/ Bcl-xL
Target	BCL2L1/Bcl-xL
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	26 kDa
Uniprot ID	Q07817
Expiration Date	12 months from date of receipt.