

Product Datasheet

Cystatin-C(8F14) antibody (orb526669)

Description	Cystatin-C(8F14) antibody
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr, IHC-P, WB
Immunogen	Recombinant human Cystatin-C
Target	Cystatin-C(8F14)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
MW	13 kDa
Uniprot ID	P01034
Expiration Date	12 months from date of receipt.