

Product Datasheet

RAD51 antibody (orb526670)

Description	RAD51 antibody
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombination human RAD51 protein
Target	RAD51
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
MW	37 kDa
Uniprot ID	Q06609
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.