


Product Datasheet

LAMP1 antibody (orb526686)

Description	LAMP1 antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from human LAMP1
Target	LAMP1
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG1,k
Clonality	Monoclonal
MW	42 kDa
Uniprot ID	P11279
Dilution Range	WB=1:500-1000
Expiration Date	12 months from date of receipt.



Western blot: Lane 1: Human HT1080 cell; La...