

Product Datasheet

Solute Carrier Family 1 antibody (orb526711)

Description	Solute Carrier Family 1 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from mouse Solute carrier family 1 (161-260/541 aa)
Target	SLC1A5
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	56 kDa
Uniprot ID	Q15758
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.