

Product Datasheet

A1-42 (C-t) antibody (orb526752)

Description	A1-42 (C-t) antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human A β 1-42 (CT) (37-42/42 aa)
Target	A β 1-42 (CT)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	4.4 kDa
Uniprot ID	P05067
Dilution Range	IHC-P=1:100-500
Expiration Date	12 months from date of receipt.