


Product Datasheet

Human CA1/Carbonic anhydrase 1 ELISA Kit (orb526784)

Description	Human CA1/Carbonic anhydrase 1 ELISA Kit
Reactivity	Human
Range	312 pg/ml - 20,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: A2-F261
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum, plasma (heparin, EDTA, citrate) and urine.
Uniprot ID	P00915
Sensitivity	50 pg/ml
Dilution Range	Dilution ratio of 1:100, concentration in serum and plasma is around 600 ng/ml.
Expiration Date	Please enquire.



Human CA1
PicoKine
ELISA
Kitstandard
cur...