


Product Datasheet

Human CA2 / Carbonic Anhydrase 2 ELISA Kit (orb526785)

Description	Human CA2 / Carbonic Anhydrase 2 ELISA Kit
Reactivity	Human
Range	0.78 ng/ml - 50 ng/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: S2-K260
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA, citrate) and urine.
Uniprot ID	P00918
Sensitivity	50 pg/ml
Dilution Range	Dilution ratio of 1:10, concentration in serum and plasma is 15-76 ng/ml.
Expiration Date	Please enquire.



Human CA2
PicoKine
ELISA
Kitstandard
cur...