


Product Datasheet

Human CD14 ELISA Kit (orb526786)

Description	Human CD14 ELISA Kit
Reactivity	Human
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: CHO; Immunogen sequence: T20-C352
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA, citrate).
Uniprot ID	P08571
Sensitivity	50 pg/ml
Dilution Range	Dilution ratio of 1:4000, concentration in serum and plasma is around 4800 ng/ml.
Expiration Date	Please enquire.



Human CD14
PicoKine
ELISA Kit
standard c...