


Product Datasheet

Human Annexin A1 ELISA Kit (orb526798)

Description	Human Annexin A1 ELISA Kit
Reactivity	Human
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: A2-N346
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA, citrate).
Uniprot ID	P04083
Sensitivity	100 pg/ml
Dilution Range	For serum, dilution ratio is 1:10 and concentration is around 20 ng/ml. For plasma, dilution ratio is 1:500 and concentration is around 1 ug/ml.
Expiration Date	Please enquire.



Human
Annexin A1
PicoKine
ELISA Kit
stan...