


Product Datasheet

Human Ki67 ELISA Kit (orb526803)

Description	Human Ki67 ELISA Kit
Reactivity	Human
Range	7.8 pg/ml - 500 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: N3120-I3256
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA, citrate).
Uniprot ID	P46013
Sensitivity	1 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is less than the lowest standard, 7.8 pg/ml.
Expiration Date	Please enquire.



Human Ki67
PicoKine
ELISA Kit
standard c...