


Product Datasheet

Human Endorepellin/HSPG2 ELISA Kit (orb526804)

Description	Human Endorepellin/HSPG2 ELISA Kit
Reactivity	Human
Range	0.624 ng/ml - 40 ng/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: E3687-S4391
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA, citrate).
Uniprot ID	P98160
Sensitivity	0.05 ng/ml
Dilution Range	Dilution ratio of 1:5, concentration in serum and plasma is 6-30 ng/ml.
Expiration Date	Please enquire.



Human
Endorepellin/HSPG2
PicoKine ELISA ...