


Product Datasheet

Human CD164 ELISA Kit (orb526811)

Description	Human CD164 ELISA Kit
Reactivity	Human
Range	93.8 pg/ml - 6,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: D24-D162
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA, citrate).
Uniprot ID	Q04900
Sensitivity	50 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is 1-2 ng/ml.
Expiration Date	Please enquire.



Human
CD164
PicoKine
ELISA Kit
standard ...