


Product Datasheet

Human REG1A ELISA Kit (orb526812)

Description	Human REG1A ELISA Kit
Reactivity	Human
Range	62.5 pg/ml - 4,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: Q23-N166
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin).
Uniprot ID	P05451
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:100, concentration in serum and plasma is around 60 ng/ml.
Expiration Date	Please enquire.



Human
REG1A
PicoKine
ELISA
Kitstandard
C...