


Product Datasheet

Rat Carboxylesterase 3/CES3/Ces1d ELISA Kit (orb526816)

Description	Rat Carboxylesterase 3/CES3/Ces1d ELISA Kit
Reactivity	Rat
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: Y19-L565
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA, citrate).
Uniprot ID	P16303
Sensitivity	50 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is around 2 ng/ml.
Expiration Date	Please enquire.



Rat ES3
PicoKine
ELISA Kit
standard
curv...