


Product Datasheet

Mouse Plexin A1/PLXNA1 ELISA Kit (orb526817)

Description	Mouse Plexin A1/PLXNA1 ELISA Kit
Reactivity	Mouse
Range	0.312 ng/ml - 20 ng/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: S28-P1242
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA).
Uniprot ID	P70206
Sensitivity	0.1 ng/ml
Dilution Range	Dilution ratio of 1:5, concentration in serum and plasma is 7.5-10 ng/ml.
Expiration Date	Please enquire.



Mouse
Plexin A1
PicoKine
ELISA Kit
stand...