


Product Datasheet

Human Plexin B3 ELISA Kit (orb526820)

Description	Human Plexin B3 ELISA Kit
Reactivity	Human
Range	93.6 pg/ml - 6,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: H45-Q1255
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA).
Uniprot ID	Q9ULL4
Sensitivity	15 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is around 2.4 ng/ml.
Expiration Date	Please enquire.



Human
Plexin B3
PicoKine
ELISA Kit
stand...