


Product Datasheet

Human REG3A ELISA Kit (orb526822)

Description	Human REG3A ELISA Kit
Reactivity	Human
Range	0.78 ng/ml - 50 ng/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: E27-D175
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin).
Uniprot ID	Q06141
Sensitivity	0.1 ng/ml
Dilution Range	Dilution ratio of 1:50, concentration in serum and plasma is 300-500 ng/ml.
Expiration Date	Please enquire.



Human
REG3A
PicoKine
ELISA Kit
standard ...