


Product Datasheet

Human LYPD3 ELISA Kit (orb526825)

Description	Human LYPD3 ELISA Kit
Reactivity	Human
Range	0.312 ng/ml - 20 ng/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: L31-C326
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA, citrate).
Uniprot ID	O95274
Sensitivity	0.1 ng/ml
Dilution Range	Dilution ratio of 1:10, concentration in serum and plasma is 5-8 ng/ml.
Expiration Date	Please enquire.



Human
LYPD3
PicoKine
ELISA Kit
standard ...