


Product Datasheet

Human CPA2 ELISA Kit (orb526832)

Description	Human CPA2 ELISA Kit
Reactivity	Human
Range	93.6 pg/ml - 6,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: S115-Y419
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA).
Uniprot ID	P48052
Sensitivity	15 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is around 3 ng/ml.
Expiration Date	Please enquire.



Human CPA2
PicoKine
ELISA Kit
standard c...