


## Product Datasheet

### Human CTRP9 ELISA Kit (orb526837)

|                            |                                                                                                                                                                                   |
|----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Human CTRP9 ELISA Kit                                                                                                                                                             |
| <b>Reactivity</b>          | Human                                                                                                                                                                             |
| <b>Range</b>               | 46.9 pg/ml - 3,000 pg/ml                                                                                                                                                          |
| <b>Tested Applications</b> | ELISA                                                                                                                                                                             |
| <b>Immunogen</b>           | Expression system for standard: NS0; Immunogen sequence: Q20-P333                                                                                                                 |
| <b>Storage</b>             | Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.) |
| <b>Note</b>                | For research use only                                                                                                                                                             |
| <b>Sample Types</b>        | cell culture supernatants, serum and plasma (heparin).                                                                                                                            |
| <b>Uniprot ID</b>          | <b>P0C862</b>                                                                                                                                                                     |
| <b>Sensitivity</b>         | 10 pg/ml                                                                                                                                                                          |
| <b>Dilution Range</b>      | For serum, dilution ratio is 1:1 and concentration is around 400 pg/ml. For plasma, dilution ratio is 1:1 and concentration is around 150 pg/ml.                                  |
| <b>Expiration Date</b>     | Please enquire.                                                                                                                                                                   |



Human  
CTRP9  
PicoKine  
ELISA Kit  
standard ...