


Product Datasheet

Mouse Plexin B2/PLXNB2 ELISA Kit (orb526838)

Description	Mouse Plexin B2/PLXNB2 ELISA Kit
Reactivity	Mouse
Range	0.78 ng/ml - 50 ng/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: CHO; Immunogen sequence: L20-P1201
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA).
Uniprot ID	B2RXS4
Sensitivity	50 pg/ml
Dilution Range	Dilution ratio of 1:10, concentration in serum and plasma is 63-100 ng/ml.
Expiration Date	Please enquire.



Mouse
Plexin B2
PicoKine
ELISA Kit
stand...