


Product Datasheet

Mouse SerpinE2 ELISA Kit (orb526840)

Description	Mouse SerpinE2 ELISA Kit
Reactivity	Mouse
Range	0.78 ng/ml - 50 ng/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: S20-P397
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA, citrate) ,urine.
Uniprot ID	Q07235
Sensitivity	0.1 ng/ml
Dilution Range	Dilution ratio of 1:10000, concentration in serum and plasma is 30-70 ug/ml.
Expiration Date	Please enquire.



Mouse
SerpinE2
PicoKine
ELISA Kit
standa...