


Product Datasheet

Rat CPE/Carboxypeptidase E ELISA Kit (orb526844)

Description	Rat CPE/Carboxypeptidase E ELISA Kit
Reactivity	Rat
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: L43-F476
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA).
Uniprot ID	P15087
Sensitivity	50 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is 800-1000 pg/ml.
Expiration Date	Please enquire.



Rat CPE
PicoKine
ELISA Kit
standard
curv...