


Product Datasheet

Human P4HB ELISA Kit (orb526845)

Description	Human P4HB ELISA Kit
Reactivity	Human
Range	93.8 pg/ml - 6,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: D18-L508
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA).
Uniprot ID	P07237
Sensitivity	15 pg/ml
Dilution Range	Dilution ratio of 1:5, concentration in serum and plasma is 2000-4000 pg/ml.
Expiration Date	Please enquire.



Human P4HB
PicoKine
ELISA Kit
standard c...