


Product Datasheet

Rat GAS6 ELISA Kit (orb526849)

Description	Rat GAS6 ELISA Kit
Reactivity	Rat
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: D115-P673
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA).
Uniprot ID	Q63772
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:5, concentration in serum and plasma is 1.2-4 ng/ml.
Expiration Date	Please enquire.



Rat GAS6
PicoKine
ELISA Kit
standard
cur...