


## Product Datasheet

### Human C1S ELISA Kit (orb526855)

|                            |                                                                                                                                                                                   |
|----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Human C1S ELISA Kit                                                                                                                                                               |
| <b>Reactivity</b>          | Human                                                                                                                                                                             |
| <b>Range</b>               | 0.78 ng/ml - 50 ng/ml                                                                                                                                                             |
| <b>Tested Applications</b> | ELISA                                                                                                                                                                             |
| <b>Immunogen</b>           | Expression system for standard: NS0; Immunogen sequence: E16-D688                                                                                                                 |
| <b>Storage</b>             | Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.) |
| <b>Note</b>                | For research use only                                                                                                                                                             |
| <b>Sample Types</b>        | cell culture supernatants, serum, plasma (heparin, EDTA or citrate) , urine, human milk                                                                                           |
| <b>Uniprot ID</b>          | <b>P09871</b>                                                                                                                                                                     |
| <b>Sensitivity</b>         | 0.1 ng/ml                                                                                                                                                                         |
| <b>Dilution Range</b>      | Dilution ratio of 1:5000, concentration in serum and plasma is 20-50 ug/ml.                                                                                                       |
| <b>Expiration Date</b>     | Please enquire.                                                                                                                                                                   |



Human C1S  
PicoKine  
ELISA  
Kitstandard  
cur...