


Product Datasheet

Human FCRL5 ELISA Kit (orb526856)

Description	Human FCRL5 ELISA Kit
Reactivity	Human
Range	0.78 ng/ml - 50 ng/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: Q16-R844
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA).
Uniprot ID	Q96RD9
Sensitivity	0.1 ng/ml
Dilution Range	Dilution ratio of 1:5, concentration in serum and plasma is around 50 ng/ml.
Expiration Date	Please enquire.



Human
FCRL5
PicoKine
ELISA
Kitstandard
C...