


Product Datasheet

Human CA4/Carbonic anhydrase 4 ELISA Kit (orb526860)

Description	Human CA4/Carbonic anhydrase 4 ELISA Kit
Reactivity	Human
Range	312 pg/ml - 20,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: A19-S284
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA, citrate).
Uniprot ID	P22748
Sensitivity	50 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is around 2 ng/ml.
Expiration Date	Please enquire.



Human CA4
PicoKine
ELISA
Kitstandard
cur...