


Product Datasheet

Human CA6/Carbonic anhydrase 6 ELISA Kit (orb526862)

Description	Human CA6/Carbonic anhydrase 6 ELISA Kit
Reactivity	Human
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: Q18-N308
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA).
Uniprot ID	P23280
Sensitivity	50 pg/ml
Dilution Range	Dilution ratio of 1:10, concentration in serum and plasma is 15-24 ng/ml.
Expiration Date	Please enquire.



Human CA6
PicoKine
ELISA
Kitstandard
cur...