


Product Datasheet

Human CD55 ELISA Kit (orb526869)

Description	Human CD55 ELISA Kit
Reactivity	Human
Range	312 pg/ml - 20,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: D35-S353
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA).
Uniprot ID	P08174
Sensitivity	50 pg/ml
Dilution Range	Dilution ratio of 1:50, concentration in serum and plasma is 20-50 ng/ml.
Expiration Date	Please enquire.



Human CD55
PicoKine
ELISA
Kitstandard
CU...