


Product Datasheet

Human Complement C1R ELISA Kit (orb526872)

Description	Human Complement C1R ELISA Kit
Reactivity	Human
Range	1.56 ng/ml - 100 ng/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: S18-D705
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA, citrate) and urine.
Uniprot ID	P00736
Sensitivity	0.1 ng/ml
Dilution Range	Dilution ratio of 1:500, concentration in serum and plasma is around 10 ug/ml.
Expiration Date	Please enquire.



Human C1R
PicoKine
ELISA
Kitstandard
cur...