


Product Datasheet

Mouse C1R ELISA Kit (orb526873)

Description	Mouse C1R ELISA Kit
Reactivity	Mouse
Range	0.625 ng/ml - 40 ng/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: S18-D705
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA, citrate) and urine.
Uniprot ID	Q8CFG9, Q8CG16
Sensitivity	0.1 ng/ml
Dilution Range	Dilution ratio of 1:10000, concentration in serum and plasma is around 100 ug/ml.
Expiration Date	Please enquire.



Mouse C1R
PicoKine
ELISA
Kitstandard
cur...