


Product Datasheet

Human CNDP1/Carnosine Dipeptidase 1 ELISA Kit (orb526875)

Description	Human CNDP1/Carnosine Dipeptidase 1 ELISA Kit
Reactivity	Human
Range	6.25 ng/ml - 400 ng/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: S27-H507
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA).
Uniprot ID	Q96KN2
Sensitivity	50 pg/ml
Dilution Range	Dilution ratio of 1:100, concentration in serum and plasma is 1.49-7.0 ug/ml.
Expiration Date	Please enquire.



Human
CNDP1
PicoKine
ELISA
Kitstandard
C...