


Product Datasheet

Mouse CD200R1 ELISA Kit (orb526876)

Description	Mouse CD200R1 ELISA Kit
Reactivity	Mouse
Range	78 pg/ml - 5,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: T26-P238
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA).
Uniprot ID	Q9ES57
Sensitivity	15 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is around 100 pg/ml.
Expiration Date	Please enquire.



Mouse
CD200R1
PicoKine
ELISA
Kitstandard...