


Product Datasheet

Human Nectin-2/CD112 ELISA Kit (orb526877)

Description	Human Nectin-2/CD112 ELISA Kit
Reactivity	Human
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: Q32-G360
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA).
Uniprot ID	Q92692
Sensitivity	0.1 ng/ml
Dilution Range	Dilution ratio of 1:10, concentration in serum and plasma is 4-10 ng/ml.
Expiration Date	Please enquire.



Human
Nectin-
2/CD112
PicoKine
ELISA Kits...